


TOPSPIN quick guide to record routine 1D ^1H , $^{13}\text{C}\{^1\text{H}\}$ and $^{13}\text{C}\{^1\text{H}\}$ DEPT135 NMR

Penn Chemistry NMR Facility

Getting Started:

- Login using your group ID and password.
 - In **Desktop**, double click on  **TOPSPIN icon**.
- BSMS keypad Operation:
- press **LOCK ON/OFF** button to switch off lock hold.
 - press **SPIN ON/OFF** button to stop spinning (if spin indicator is **ON**).
 - press **LIFT ON/OFF** button to eject sample.
- Gently slide NMR tube into the spinner and adjust tube height in depth gauge before inserting sample.
- On BSMS keypad, press **LIFT ON/OFF** button to insert sample.

Creating New Dataset:

- Type **new** ↵
- Enter **sample_id** (maximum 13 character) in **NAME** box, type a numeric value for **EXPNO** (1-999) and a numeric value for **PROCNO** (1-999).
- Click the down-arrow of the **Solvent** box to select a solvent
- Click the down-arrow of the **Experiment** box to select **use current params**
- Type **rpar** ↵ and select parameter set
- Parameter Name: proton = **1_protonstd**
carbon = **1_carbonstd** DEPT135 = **1_dept135**

Locking and Shimming:



- Type **lopo** ↵. Select deuterated solvent from the list.
- Type **lockdisp** ↵.
- If lock ring down pattern is normal, press **LOCK ON/OFF** button to start lock hold.
- Optimize **X** and **Y** shim iteratively without spinning the sample.
- Press **SPIN ON/OFF** button to start sample spinning (wait for spin indicator to stop blinking).
- Optimize **Z** and **Z2** shim iteratively (make sure **FINE** indicator is **ON**)
- Press **STDBY**.

Acquiring Spectrum:


- Type **rga** ↵ and wait for **rga:finished** message.
- Type **zg** ↵ to start acquisition.
- After 4-8 transients (**only for ^1H**), type **tr** ↵ to transfer data to host computer.
- Type **efp** ↵ to perform Fourier transformation including apodization and default phasing.
- Type **apk** ↵ to perform automatic phase correction.
- If signal-to-noise of the spectrum is satisfactory, then type **halt** ↵, wait for **zg:acquisition finished** message.
- Type **efp** ↵ to transform rest of the acquired transients.

Processing Spectrum:



Phase correction:

If phase of the spectrum is fine, go to next step, otherwise click on the  button, click-hold the **0** and **1** buttons to get correct phase. Click 




Chemical shift calibration:

Horizontally expand solvent peak region and click the  button, click on top of the solvent peak, enter **correct ppm value** of the solvent from the chart.

Integration:

Click on the  button, with the green highlighted  button, select integration region by dragging the mouse to cover the whole region of the peak(s) [**Adjust bias and slope if necessary**]. Select the known integral by right-clicking and choosing **Calibrate** from the popup menu while cursor is under the integral region, enter value for the integral in the popup box. Click **OK** button.

Peak picking:

Click on the  button, with the green highlighted  button, select peak picking range, and with the  button, modify the peak picking range.

Printing Spectrum:

- Click **Title** tab and enter plot title
- Type **layout** ↵ and select **desired layout** by clicking down-arrow button in **LAYOUT** box
- Type **plot** ↵, if necessary, modify layout in TOPSPIN plot editor
- Click **File**→**print** to print the spectrum in TOPSPIN plot editor.

Logout Steps:

- On BSMS keypad:
 - press **LOCK ON/OFF** button to switch off lock hold.
 - press **SPIN ON/OFF** button to stop spinning.
 - press **LIFT ON/OFF** button to eject sample.
- Insert standard sample tube back into the magnet.
- Start Lock hold by typing **lopo cdcl3** ↵ (**don't spin standard sample**)
- Optimize **Z1** and **Z2** shims without spinning.
- Click **quit** tab in **Lockdisp** window to close.
- Type **exit** ↵ to close **TOPSPIN** window.
- In **desktop**, Click on **start**, select **Log Off**, click **Log Off** tab to log out.
- Fill out log book with all necessary info (your ID, group ID, solvent, time, experiment).

Emergency contact: Dr. George Furst, furst@sas.upenn.edu, 215-898-3407
Dr. Jun Gu, jungu@sas.upenn.edu, 215-898-1302